



1. Introduction

Olink® IMMUNO-ONCOLOGY is a reagent kit that measures 92 immuno-oncology related human proteins simultaneously using just 1 µL of serum, plasma or other human sample type. The analytical performance of the product has been carefully validated and the results are presented in this document. Please note that when a new panel is developed, both the individual assays and 92-plex panel as a whole are subject to our thorough validation procedure. If individual assays are subsequently improved or one or more assays are replaced in later versions of the panel, focus is placed on thoroughly validating the individual assays in question.

1.1 TECHNOLOGY

The Olink reagents are based on the Proximity Extension Assay (PEA) technology^{1,2}, where 92 oligonucleotide-labeled antibody probe pairs are allowed to bind to their designated target protein, if it is present in the sample. A PCR reporter sequence is formed by a proximity-dependent DNA hybridization and polymerization event. This is then measured, using real-time PCR. The assay is performed in a homogeneous 96-well format without any need for washing steps, see Figure 1.

1.2 QUALITY CONTROLS

Internal and external controls have been developed by Olink for data normalization and quality control. They have been designed to monitor the technical performance of the assay, as well as the quality of individual samples. This provides information at each step of the Olink protocol (see Figure 1). The internal controls are added to each sample and include two Incubation controls, one Extension control and one Detection control. The Incubation controls (two non-human proteins) monitor all three steps starting with the immuno reaction. The Extension Control (an antibody linked to two matched oligonucleotides for immediate proximity-dependent hybridization and extension that does not require

antibody binding to the target protein) monitors the extension and detection steps and is used for data normalization across samples. Finally, the Detection control (a synthetic double-stranded template) specifically monitors the detection step. If one or more of the internal control values deviate from a pre-determined range, the sample will be flagged and may be removed before statistical analysis. An external control called the inter-plate control (IPC), is included on each plate and used in a second normalization step. The IPC is made up of a pool of probes similar to the Extension control (Ext Ctrl), except that it is generated with 92 matching oligonucleotide pairs. Furthermore, the IPC improves inter-assay precision and allows for optimal comparison of data derived from multiple runs. The term “Normalized Protein eXpression (NPX)” refers to data that has been normalized as described above.

1.3 DATA ANALYSIS

The data analysis described in this document was performed using a pre-processing normalization procedure. For each sample and data point, the corresponding Cq-value for the Extension control was subtracted, normalizing for technical variation within one run. Normalization between runs was then performed for each assay by subtracting the corresponding dCq-value for the Interplate Control (IPC) from the dCq-values generated. In the final step of the pre-processing procedure the values were set relative to a correction factor that is determined by Olink for each batch of conjugated PEA probes. The Normalized Protein eXpression (NPX) unit generated by these procedures is on a log₂ scale where a higher number represents a higher level of the target in the sample, typically with the background level at around zero. Linearization of data is performed by the mathematical operation 2^{NPX}. Coefficient of variation (CV) calculations were performed on linearized values.




Fig 1. Olink assay procedure (above) and controls (below). The internal controls enable monitoring of the three core steps in the Olink assay and used for quality control and data normalization. Detection is performed by using the Fluidigm® Biomark™ or the Fluidigm® Biomark™ HD system.

2. Performance characteristics

2.1 SAMPLE TYPES

The ability to use different sample types with Olink IMMUNO-ONCOLOGY was evaluated by collecting matched EDTA, acid citrate dextrose (ACD), and sodium heparin plasma and serum samples from 4 healthy individuals. Table 1 summarizes response values for 32 normal EDTA plasma samples expressed in NPX, as well as the relative differences for the other sample types compared to EDTA plasma. Variations observed between responses in the different sample types tested were generally small, and all assays will therefore function without limitation in these sample types. In addition, cell lysates from 10 different cell lines were also evaluated. For more information about using cell lysates please contact support@olink.com.

2.2 ANALYTICAL MEASUREMENT

NOTE: The technical performance data based on *in vitro* assays using recombinant antigen must **NOT** be used to derive actual concentrations of native proteins in biological samples from the relative quantification NPX data that is obtained from an Olink assay.

DETECTION LIMIT

Calibrator curves were determined for all biomarkers, for which recombinant antigen was available, simultaneously in a multiplex format. For these assays, the Limit of detection (LOD) was defined as 3 standard deviations above background and reported in pg/mL (see Table 1 and Figure 2).

HIGH DOSE HOOK EFFECT

The high dose hook effect is a state of antigen excess relative to the reagent antibodies, resulting in falsely lower values reported, which can lead to misinterpretation of results. Therefore, the hook effect was determined for each analyte where applicable, and reported in pg/mL (see Table 1).

MEASURING RANGE

The analytical measuring range was defined by the lower limit of quantification (LLOQ) and upper limit of quantification (ULOQ) and reported in log₁₀, see Table 1. The upper and lower limits of quantification, ULOQ and LLOQ, respectively were calculated with the following trueness and precision criteria; relative error ≤ 30% and CV ≤ 30% of back-calculated values, and reported in pg/mL (see Table 1).

Three selected assays with their analytical data are shown in Figure 2 and the distribution of measuring ranges of 90 assays compared to endogenous plasma levels are shown in Figure 3. Where applicable, individual calibrator curves are available on the specific biomarker page on the Olink website.




Fig 2. Calibrator curves from 3 assays and their corresponding analytical measurement data.




Fig 3. Distribution of analytical measuring range (light blue bars), defined by the lower and upper limits of quantification (LLOQ-ULOQ), and normal plasma levels (dark blue bars) for all assays with currently available data.

Table 1. Assay performance parameters. Sample Types; Normalized Protein eXpression (NPX), Endogenous Interference, Analytical rNge; Limit of Detection (LOD), Lower Limit of Quantification (LLOQ), Upper Limit of Quantification (ULOQ), High Dose Effect (Hook), Range and Precision indicative of assay performance are shown for the analytes where applicable. Not available, NA

| Target | UniProt No | Sample types | | | Endogenous interference | | | Analytical range | | | | Precision | | | |
|--|---------------|----------------------------|--------|------------|-----------------------------|---------|-------|------------------|--------|--------|----------|-----------|-------|-------|-------|
| | | Normal plasma levels (NPX) | | | Relative to EDTA plasma (%) | | | (mg/mL) | pg/mL | | | | log10 | | % CV |
| | | 10th %tile | Median | 90th %tile | ACD | Heparin | Serum | Hemolysate | LOD | LLOQ | ULOQ | Hook | Range | Intra | Inter |
| Adenosine deaminase (ADA) | P00813 | 3.3 | 3.8 | 4.5 | 94 | 90 | 103 | 0.2 | 244.1 | 244.1 | 1000000 | 1000000 | 3.6 | 5.8 | 12.7 |
| Adhesion G-protein coupled receptor G1 (ADGRG1) | Q9Y653 | NA | 1.5 | 1.9 | 94 | 86 | 120 | 15 | 1953.1 | 1953.1 | 500000 | 1000000 | 2.4 | 8.0 | 14.5 |
| Angiotensin-1 (ANGPT1) | Q15389 | 5.9 | 7.4 | 10.4 | 16 | 157 | 274 | 15 | 122.1 | 244.1 | 500000 | 500000 | 3.3 | 7.4 | 8.2 |
| Angiotensin-1 receptor (TIE2) | Q02763 | 7.6 | 8.1 | 8.5 | 97 | 92 | 105 | 15 | 244.1 | 244.1 | 1000000 | 1000000 | 3.6 | 7.2 | 8.4 |
| Angiotensin-2 (ANGPT2) | Q15123 | 4.7 | 5.3 | 6.1 | 90 | 84 | 107 | 15.0 | 488.3 | 488.3 | 1000000 | 1000000 | 3.3 | 7.4 | 12.4 |
| Arginase-1 (ARG1) | P05089 | NA | 1.6 | 4.6 | 33 | 52 | 163 | 0.0 | 244.1 | 244.1 | 31250 | 250000 | 2.1 | 6.5 | 7.0 |
| Carbonic anhydrase 9 (CAIX) | Q16790 | 3.3 | 4.1 | 4.8 | 90 | 100 | 110 | 15 | 15.3 | 30.5 | 62500 | 125000 | 3.3 | 7.7 | 7.3 |
| Caspase-8 (CASP-8) | Q14790 | 4.2 | 4.8 | 5.7 | 72 | 130 | 174 | 0.2 | 3.8 | 3.8 | 62500 | 62500 | 4.2 | 7.2 | 9.7 |
| C-C motif chemokine 13 (MCP-4) | Q99616 | 6.1 | 7.2 | 8.8 | 57 | 120 | 227 | 15 | 7.6 | 30.5 | 7812 | 7812 | 2.4 | 8.3 | 13.1 |
| C-C motif chemokine 17 (CCL17) | Q92583 | 5.7 | 7.4 | 9.4 | 37 | 134 | 276 | 15 | 1.9 | 1.9 | 3906 | 7812 | 3.3 | 8.7 | 13.9 |
| C-C motif chemokine 19 (CCL19) | Q99731 | 7.9 | 8.6 | 10.3 | 98 | 87 | 108 | 15 | 122.1 | 244.1 | 62500 | 125000 | 2.4 | 7.7 | 11.4 |
| C-C motif chemokine 2 (MCP-1) | P13500 | 9.2 | 9.8 | 10.6 | 100 | 99 | 122 | 15.0 | 0.5 | 0.5 | 3906 | 3906 | 3.9 | 8.0 | 11.9 |
| C-C motif chemokine 20 (CCL20) | P78556 | 5.1 | 6.2 | 8.2 | 102 | 75 | 74 | 15 | 15.3 | 30.5 | 15625 | 31250 | 2.7 | 8.5 | 16.0 |
| C-C motif chemokine 23 (CCL23) | P55773 | 8.7 | 9.4 | 10.6 | 95 | 85 | 90 | 15 | 30.5 | 61.0 | 31250 | 62500 | 2.7 | 8.2 | 14.1 |
| C-C motif chemokine 3 (CCL3) | P10147 | 4.6 | 5.2 | 6.4 | 70 | 98 | 145 | 15 | 0.5 | 0.5 | 7812 | 3906 | 4.2 | 8.2 | 14.5 |
| C-C motif chemokine 4 (CCL4) | P13236 | 6.3 | 7.1 | 8.3 | 66 | 99 | 134 | 15.0 | 7.6 | 7.6 | 62500.0 | 125000 | 3.9 | 8.5 | 15.7 |
| C-C motif chemokine 7 (MCP-3) | P80098 | 1.6 | 2.3 | 3.0 | 85 | 111 | 127 | 7.5 | 30.5 | 61.0 | 125000.0 | 1000000 | 3.3 | 8.7 | 10.2 |
| C-C motif chemokine 8 (MCP-2) | P80075 | 5.4 | 6.9 | 8.3 | 62 | 92 | 173 | 15 | 1.9 | 3.8 | 3906 | 3906 | 3.0 | 7.6 | 12.7 |
| CD27 antigen (CD27) | P26842 | 7.6 | 8.2 | 8.7 | 92 | 100 | 108 | 15 | 7.6 | 15.3 | 62500 | 125000 | 3.6 | 6.7 | 7.2 |
| CD40 ligand (CD40-L) | P29965 | 2.9 | 4.3 | 7.0 | 39 | 292 | 1042 | 15 | 7.6 | 7.6 | 15625 | 31250 | 3.3 | 7.7 | 9.4 |
| CD40L receptor (CD40) | P25942 | 10.2 | 10.7 | 11.2 | 95 | 105 | 126 | 15.0 | 0.1 | 0.1 | 7812 | 31250 | 4.8 | 7.3 | 10.0 |
| CD70 antigen (CD70) | P32970 | 3.4 | 4.0 | 4.7 | 92 | 97 | 169 | 15 | 122.1 | 122.1 | 62500 | 125000 | 2.7 | 8.0 | 11.2 |
| CD83 antigen (CD83) | Q01151 | 2.5 | 3.0 | 3.5 | 88 | 83 | 104 | 15 | 3.8 | 15.3 | 31250 | 62500 | 3.3 | 7.8 | 11.0 |
| C-X-C motif chemokine 1 (CXCL1) | P09341 | 7.2 | 8.7 | 10.8 | 37 | 157 | 285 | 7.5 | 7.6 | 30.5 | 31250 | 31250 | 3.0 | 7.2 | 11.8 |
| C-X-C motif chemokine 10 (CXCL10) | P02778 | 6.7 | 7.5 | 9.5 | 88 | 87 | 112 | 15.0 | 30.5 | 61.0 | 31250 | 31250 | 2.7 | 11.5 | 14.3 |
| C-X-C motif chemokine 11 (CXCL11) | Q14625 | 4.6 | 5.8 | 8.9 | 29 | 147 | 281 | 0.5 | 30.5 | 61.0 | 31250 | 31250 | 2.7 | 8.8 | 12.2 |
| C-X-C motif chemokine 13 (CXCL13) | Q43927 | 8.1 | 8.8 | 9.8 | 106 | 81 | 131 | 15 | 30.5 | 30.5 | 15625 | 15625 | 2.7 | 6.2 | 10.6 |
| C-X-C motif chemokine 5 (CXCL5) | P42830 | 7.8 | 10.2 | 13.1 | 13 | 190 | 271 | 7.5 | 7.6 | 7.6 | 15625 | 31250 | 3.3 | 7.9 | 11.5 |
| C-X-C motif chemokine 9 (CXCL9) | Q07325 | 6.2 | 6.9 | 8.2 | 89 | 94 | 95 | 15 | 1.9 | 7.6 | 31250 | 31250 | 3.6 | 9.8 | 13.6 |
| Cytotoxic and regulatory T-cell molecule (CRTAM) | Q95727 | 4.1 | 4.8 | 5.8 | 90 | 105 | 127 | 15.0 | 30.5 | 61.0 | 125000 | 125000 | 3.3 | 7.2 | 15.1 |
| Decorin (DCN) | P07585 | 4.7 | 5.0 | 5.4 | 99 | 120 | 118 | 15.0 | 488.3 | 488.3 | 31250 | 1000000 | 1.8 | 8.2 | 8.2 |
| Fibroblast growth factor 2 (FGF2) | P09038 | NA | NA | 2.0 | 95 | 56 | 89 | 15 | 30.5 | 122.1 | 250000 | 1000000 | 3.3 | 6.3 | 9.0 |
| Fractalkine (CX3CL1) | P78423 | 6.1 | 6.6 | 7.1 | 94 | 110 | 143 | 15 | 30.5 | 61.0 | 15625 | 1000000 | 2.4 | 7.1 | 7.5 |
| Galectin-1 (Gal-1) | P09382 | 6.4 | 6.7 | 7.2 | 94 | 102 | 111 | 15 | 244.1 | 244.1 | 125000.0 | 250000 | 2.7 | 6.0 | 8.0 |
| Galectin-9 (Gal-9) | Q00182 | 7.7 | 8.1 | 8.5 | 102 | 102 | 110 | 0.5 | 3.8 | 7.6 | 3906 | 15625 | 2.7 | 5.3 | 9.1 |
| Granzyme A (GZMA) | P12544 | 5.1 | 5.5 | 6.2 | 90 | 88 | 109 | 7.5 | 30.5 | 30.5 | 62500 | 125000 | 3.3 | 9.9 | 13.4 |
| Granzyme B (GZMB) | P10144 | 2.2 | 3.0 | 3.6 | 89 | 76 | 88 | 0.5 | 1.9 | 1.9 | 15625 | 31250 | 3.9 | 8.1 | 13.4 |
| Granzyme H (GZMH) | P20718 | 3.5 | 4.5 | 5.7 | 104 | 103 | 116 | 0.5 | 30.5 | 61.0 | 62500 | 125000 | 3.0 | 8.5 | 12.1 |
| Heme oxygenase 1 (HO-1) | P09601 | 11.5 | 12.2 | 12.8 | 94 | 92 | 98 | 15.0 | 15.3 | 15.3 | 31250 | 31250 | 3.3 | 7.3 | 13.1 |
| Hepatocyte growth factor (HGF) | P14210 | 6.9 | 7.5 | 8.3 | 74 | 73 | 153 | 15 | 3.8 | 15.3 | 31250 | 250000 | 3.3 | 8.4 | 10.0 |
| ICOS ligand (ICOSLG) | Q75144 | 5.3 | 5.7 | 6.0 | 95 | 144 | 146 | 15 | 61.0 | 61.0 | 62500 | 1000000 | 3.0 | 6.6 | 7.4 |
| Interferon gamma (IFN-gamma) | P01579 | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* |
| Interleukin-1 alpha (IL-1 alpha) | P01583 | NA | NA | NA | 88 | 81 | 69 | 15.0 | NA | NA | NA | NA | NA | 19.6 | NA |
| Interleukin-10 (IL10) | P22301 | 2.1 | 2.6 | 3.5 | 87 | 85 | 111 | 15 | 15.3 | 30.5 | 1000000 | 1000000 | 4.5 | 9.7 | 12.0 |
| Interleukin-12 (IL12) | P29459,P29460 | 4.9 | 5.8 | 6.8 | 84 | 75 | 109 | 15.0 | 3.8 | 3.8 | 31250 | 62500 | 3.9 | 7.5 | 11.4 |
| Interleukin-12 receptor subunit beta-1 (IL12RB1) | P42701 | 2.0 | 2.3 | 2.8 | 91 | 84 | 98 | 15 | 30.5 | 30.5 | 125000 | 125000 | 3.6 | 7.1 | 11.1 |
| Interleukin-13 (IL13) | P35225 | NA | NA | 0.4 | 93 | 64 | 110 | 15 | 244.1 | 488.3 | 125000 | 1000000 | 2.4 | 8.9 | 6.6 |
| Interleukin-15 (IL15) | P40933 | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* |

| Target | UniProt No | Sample types | | | | | | Endogenous interference | Analytical range | | | | Precision | | |
|---|---------------|----------------------------|--------|------------|-----------------------------|---------|-------|-------------------------|------------------|--------|---------|---------|-----------|-------|-------|
| | | Normal plasma levels (NPX) | | | Relative to EDTA plasma (%) | | | (mg/mL) | pg/mL | | | log10 | % CV | | |
| | | 10th %tile | Median | 90th %tile | ACD | Heparin | Serum | Hemolysate | LOD | LLOQ | ULOQ | Hook | Range | Intra | Inter |
| Interleukin-18 (IL18) | Q14116 | 8.0 | 8.9 | 9.6 | 87 | 92 | 108 | 3.8 | 1.0 | 3.8 | 62500 | 62500 | 4.2 | 8.0 | 12.5 |
| Interleukin-2 (IL2) | P60568 | NA | NA | NA | 107 | 60 | 122 | 15 | 122.1 | 122.1 | 500000 | 1000000 | 3.6 | 11.8 | 7.2 |
| Interleukin-33 (IL33) | O95760 | NA | NA | 0.8 | 102 | 77 | 112 | 15 | 244.1 | 244.1 | 62500 | 125000 | 2.4 | 7.3 | 9.8 |
| Interleukin-4 (IL4) | P05112 | NA | NA | NA | 100 | 76 | 105 | 15 | 7.6 | 7.6 | 62500 | 62500 | 3.9 | 7.1 | 13.3 |
| Interleukin-5 (IL5) | P05113 | NA | 0.9 | 5.2 | 94 | 94 | 102 | 15.0 | 7.6 | 7.6 | 62500 | 125000 | 3.9 | 11.9 | 15.4 |
| Interleukin-6 (IL6) | P05231 | 3.4 | 4.1 | 5.2 | 106 | 106 | 119 | 15.0 | 1.0 | 1.9 | 31250 | 31250 | 4.2 | 7.3 | 10.8 |
| Interleukin-7 (IL7) | P13232 | 3.3 | 4.2 | 6.1 | 45 | 102 | 307 | 15.0 | 1.9 | 3.8 | 15625 | 31250 | 3.6 | 7.6 | 10.9 |
| Interleukin-8 (IL8) | P10145 | 4.5 | 5.3 | 7.1 | 71 | 116 | 175 | 7.5 | 0.5 | 1.0 | 3906 | 7812 | 3.6 | 7.9 | 13.5 |
| Killer cell immunoglobulin-like receptor 3DL1 (KIR3DL1) | P43629 | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* |
| Latency-associated peptide transforming growth factor beta-1 (LAP TGF-beta-1) | P01137 | 1.6 | 1.9 | 3.1 | 77 | 102 | 168 | 15.0 | 1953.1 | 1953.1 | 250000 | 250000 | 2.1 | 6.4 | 13.3 |
| Lymphocyte activation gene 3 protein (LAG3) | P18627 | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* |
| Lysosome-associated membrane glycoprotein 3 (LAMP3) | Q9UQV4 | 4.2 | 5.9 | 6.7 | 90 | 94 | 111 | 15 | 30.5 | 61.0 | 62500 | 250000 | 3.0 | 7.2 | 9.6 |
| Macrophage colony-stimulating factor 1 (CSF-1) | P09603 | 7.9 | 8.2 | 8.5 | 94 | 92 | 111 | 15 | 0.2 | 0.5 | 7812 | 31250 | 4.2 | 7.0 | 11.3 |
| Macrophage metalloproteinase-12 (MMP12) | P39900 | 6.1 | 7.1 | 8.0 | 138 | 117 | 125 | 15 | 30.5 | 122.1 | 250000 | 500000 | 3.3 | 8.1 | 9.8 |
| Matrix metalloproteinase-7 (MMP7) | P09237 | 6.7 | 8.4 | 9.8 | 442 | 421 | 472 | 15.0 | 7.6 | 7.6 | 31250 | 62500 | 3.6 | 7.0 | 11.3 |
| MHC class I polypeptide-related sequence A/B (MIC-A/B) | Q29983,Q29980 | NA | 4.6 | 5.5 | 88 | 95 | 109 | 15.0 | 30.5 | 30.5 | 15625 | 1000000 | 2.7 | 6.5 | 10.0 |
| Mucin-16 (MUC-16) | Q8WXI7 | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* |
| Natural cytotoxicity triggering receptor (NCR1) | O76036 | 3.7 | 4.2 | 4.9 | 90 | 97 | 118 | 15 | 7.6 | 7.6 | 31250 | 62500 | 3.6 | 7.3 | 8.9 |
| Natural killer cell receptor 2B4 (CD244) | Q9BZW8 | 5.7 | 6.1 | 6.6 | 86 | 97 | 113 | 15.0 | 1.9 | 3.8 | 31250 | 31250 | 3.9 | 7.4 | 9.8 |
| Natural killer cells antigen CD94 (KLRD1) | Q13241 | 4.7 | 5.6 | 6.8 | 84 | 94 | 109 | 15 | 3.8 | 3.8 | 31250 | 31250 | 3.9 | 6.9 | 10.0 |
| Nitric oxide synthase, endothelial (NOS3) | P29474 | NA | 0.9 | 2.4 | 128 | 59 | 76 | 15.0 | 976.6 | 1953.1 | 1000000 | 1000000 | 2.7 | 18.7 | 19.5 |
| Placenta growth factor (PGF) | P49763 | 7.9 | 8.3 | 8.8 | 88 | 95 | 109 | 15.0 | 1.0 | 1.0 | 7812 | 31250 | 3.9 | 7.7 | 11.3 |
| Platelet-derived growth factor subunit B (PDGF subunit B) | P01127 | 7.6 | 9.1 | 11.2 | 22 | 102 | 179 | 15 | 15.3 | 30.5 | 31250 | 62500 | 3.0 | 8.2 | 13.3 |
| Pleiotrophin (PTN) | P21246 | NA | 1.6 | 2.7 | 73 | 22 | 42 | 15 | 488.3 | 976.6 | 62500 | 125000 | 1.8 | 8.4 | 17.6 |
| Pro-epidermal growth factor (EGF) | P01133 | 3.9 | 5.4 | 8.4 | 28 | 145 | 710 | 15 | 1.0 | 1.0 | 3906 | 3906 | 3.6 | 8.0 | 9.6 |
| Programmed cell death 1 ligand 1 (PD-L1) | Q9NZQ7 | 4.3 | 4.9 | 5.6 | 69 | 91 | 108 | 15 | 7.6 | 7.6 | 250000 | 1000000 | 4.5 | 8.9 | 10.6 |
| Programmed cell death 1 ligand 2 (PD-L2) | Q9BQ51 | 2.2 | 2.6 | 3.0 | 93 | 92 | 111 | 15 | 244.1 | 244.1 | 250000 | 500000 | 3.0 | 6.4 | 11.1 |
| Programmed cell death protein 1 (PDCD1) | Q15116 | 3.3 | 4.0 | 4.7 | 89 | 98 | 111 | 15 | 7.6 | 7.6 | 125000 | 125000 | 4.2 | 9.7 | 12.9 |
| Stromal cell-derived factor 1 (CXCL12) | P48061 | NA | 0.6 | 1.8 | 88 | 60 | 76 | 15.0 | 7812.5 | 7812.5 | 62500 | 125000 | 0.9 | 7.6 | 13.1 |
| T-cell surface glycoprotein CD4 (CD4) | P01730 | NA | NA | NA | 109 | 88 | 81 | 15.0 | 244.1 | 244.1 | 1000000 | 1000000 | 3.6 | 7.1 | 9.6 |
| T-cell surface glycoprotein CD5 (CD5) | P06127 | 4.6 | 5.3 | 5.9 | 93 | 94 | 105 | 15 | 1.9 | 1.9 | 31250 | 31250 | 4.2 | 6.9 | 16.9 |
| T-cell surface glycoprotein CD8 alpha chain (CD8A) | P01732 | 8.4 | 9.6 | 10.5 | 100 | 86 | 89 | 15.0 | NA | NA | NA | NA | NA | 10.5 | 10.3 |
| T-cell-specific surface glycoprotein CD28 (CD28) | P10747 | 1.2 | 1.5 | 2.0 | 97 | 81 | 106 | 15 | 61.0 | 122.1 | 62500 | 125000 | 2.7 | 5.8 | 12.9 |
| TNF-related apoptosis-inducing ligand (TRAIL) | P50591 | 7.7 | 8.3 | 8.8 | 100 | 99 | 114 | 15.0 | 7.6 | 15.3 | 31250 | 31250 | 3.3 | 7.6 | 8.9 |
| Tumor necrosis factor (TNF) | P01375 | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* | NA* |
| Tumor necrosis factor ligand superfamily member 12 (TWEAK) | O43508 | 7.6 | 8.3 | 9.0 | 81 | 92 | 129 | 15.0 | 30.5 | 30.5 | 1000000 | 1000000 | 4.5 | 12.5 | 11.0 |
| Tumor necrosis factor ligand superfamily member 14 (TNFSF14) | O43557 | 3.1 | 3.8 | 4.6 | 71 | 125 | 341 | 7.5 | 15.3 | 15.3 | 31250 | 31250 | 3.3 | 7.8 | 8.9 |
| Tumor necrosis factor ligand superfamily member 6 (FASLG) | P48023 | 6.1 | 6.7 | 7.7 | 109 | 100 | 109 | 15.0 | 0.5 | 0.5 | 31250 | 31250 | 4.8 | 7.7 | 11.2 |
| Tumor necrosis factor receptor superfamily member 12A (TNFRSF12A) | Q9NP84 | 5.8 | 6.5 | 7.3 | 96 | 98 | 92 | 15 | 244.1 | 244.1 | 1000000 | 1000000 | 3.6 | 9.2 | 13.1 |
| Tumor necrosis factor receptor superfamily member 21 (TNFRSF21) | O75509 | 7.7 | 8.1 | 8.6 | 85 | 90 | 115 | 15 | 30.5 | 122.1 | 62500 | 125000 | 2.7 | 7.4 | 10.7 |
| Tumor necrosis factor receptor superfamily member 4 (TNFRSF4) | P43489 | 2.9 | 3.4 | 4.0 | 90 | 97 | 124 | 15.0 | 3.8 | 7.6 | 31250 | 31250 | 3.6 | 7.1 | 12.2 |
| Tumor necrosis factor receptor superfamily member 9 (TNFRSF9) | O07011 | 5.4 | 5.9 | 6.6 | 96 | 96 | 110 | 15 | 1.0 | 1.0 | 31250 | 31250 | 4.5 | 7.5 | 10.3 |
| Vascular endothelial growth factor A (VEGFA) | P15692 | 7.9 | 8.4 | 9.2 | 75 | 92 | 142 | 15.0 | 3.8 | 3.8 | 15625 | 31250 | 3.6 | 7.8 | 10.1 |
| Vascular endothelial growth factor receptor 2 (VEGFR-2) | P35968 | 7.1 | 7.7 | 8.1 | 100 | 96 | 111 | 15.0 | 7.6 | 7.6 | 31250 | 31250 | 3.6 | 6.9 | 12.2 |

*These recently updated assays were subject to rigorous validation and QC during their development, but final validation data in full-panel context is not yet available. This will be updated as soon as possible.

2.3 PRECISION

REPEATABILITY

Intra-assay variation (within-run) was calculated as the mean %CV for 6 different samples, run in triplicate, in 8 separate runs during the validation studies. Inter-assay variation (between runs) was calculated between experiments with the same operator. The reported inter-assay %CV is the average %CV for three different operators. These calculations were performed on linearized values for all analytes where response levels could be measured in serum and normal plasma, see Table 1.

Across all 92 assays, the mean intra-assay and inter-assay variations were observed to be 8.3% and 11.5%, respectively. The distribution of both intra-assay and inter-assay variations are shown in Figure 4.




Fig 4. Distribution of intra-assay and inter-assay variations of Olink IMMUNO-ONCOLOGY

REPRODUCIBILITY

Variations due to different operators in different laboratories using different equipment are another potential source of assay variation. Olink has Analysis Service labs in Sweden and the USA, and in addition there are many Olink-certified core laboratories around the world running the Olink platform (see www.olink.com/service for details). Our experience over several years is that inter-site reproducibility is very good providing that operators are properly trained. For more information please contact support@olink.com.

2.4 ANALYTICAL SPECIFICITY

ASSAY SPECIFICITY

The antibodies selected for use in Olink IMMUNO-ONCOLOGY have previously been evaluated against the 92 panel-specific proteins as well as against an additional 107 proteins. In principle, the specificity is tested by creating a test sample consisting of a pool of antigens, which is then incubated with all 92 antibody probe pairs from the panel. Only if there is a correct match will a reporter sequence be created and serve as a template for subsequent real-time qPCR. Ten sub-pools of antigen are evaluated to cover the 92 assays, see Figure 5.




Fig 5. Assay readout specificity of the Olink platform. For each assay, specificity is confirmed by testing antigen sub-pools against the complete 92-plex pool.

ENDOGENOUS INTERFERENCE

Endogenous interference from heterophilic antibodies, e.g. human anti-mouse antibody (HAMA), and rheumatoid factor is known to cause problems in some immunoassays. Evaluation of the potential impact of this specific interference was investigated during the validation of previous panels. No interference due to HAMA or RF was detected for any of the samples in previously tested panels, indicating sufficient blocking of these agents (data not shown).




Fig 6. Endogenous interference. Levels tested for hemolysate were 0.23-15 g/L hemoglobin. The highest hemolysate concentration translates to about 10% hemolysis.

Bilirubin, lipids and hemolysate, are plasma and serum components that are known to interfere with some analytical assays. These were evaluated for potential impact on the Olink assays at different added concentrations. An example of the hemolysate levels tested is shown in Figure 6. These additions represent different patient health conditions and/or sample collection irregularities. In 14 out of 92 assays, altered signal was observed by the addition of hemolysate. The reason is most likely due to the specific analytes leaking out of the disrupted blood cells. A concentration of 15 g/L of hemolysate represents 10% hemolysis of a sample. Table 1 reports the highest concentration of hemolysate that does not have an impact on assay performance.

Interference by bilirubin and lipids has previously been evaluated, and disturbance was only observed at extreme levels corresponding to 8 or 10 times normal^{3,4} values. This test was not therefore repeated for Olink IMMUNO-ONCOLOGY.

2.5 SCALABILITY

Assay performance has been previously evaluated with regard to scalability, meaning the capability of the Olink technology to maintain the same quality of performance irrespective of multiplex level. Previously, we have shown that a step-wise increase of multiplex grade (8, 24, 48, 72 and 96) does not compromise assay performance (data not shown). To further strengthen that Olink provides consistent results, single-plex assays for Growth Hormone (GH) and Matrix Metalloproteinase (MMP-7) were compared when run in a full 96-plex reaction. The results for each assay and their observed dCq-values were plotted against the entire 96-plex reaction. The square of the correlation coefficient (R^2) value was generated by linear regression.




Fig 7. Scalability of the Olink technology platform. The experiment was performed using the Olink CARDIOVASCULAR II panel. Human plasma samples were analyzed in singleplex for Growth Hormone (GH) and Matrix Metalloproteinase (MMP-7) with the equivalent assays performed in a full 96-plex reaction. The observed dCq (log2) values were plotted, and the correlation coefficient R^2 value was generated by linear regression.

3. References

1. Assarsson E, Lundberg M, Holmquist G, Björkesten J, Bucht Thorsen S, Ekman D, Eriksson A, Rennel Dickens E, Ohlsson S, Edfeldt G, Andersson AC, Lindstedt P, Stenvang J, Gullberg M, Fredriksson S. Homogenous 96-Plex PEA Immunoassay Exhibiting High Sensitivity, Specificity, and Excellent Scalability. *PLoS One* April (2014). doi: 10.1371/journal.pone.0095192.
2. Lundberg M, Eriksson A, Tran B, Assarsson E, Fredriksson S. Homogeneous antibody-based proximity extension assays provide sensitive and specific detection of low abundant proteins in human blood. *Nucleic Acid Res* June (2011). doi: 10.1093/nar/gkr424.
3. <http://emedicine.medscape.com/article/2074115-overview>
4. <http://www.nlm.nih.gov/medlineplus/ency/article/003479.htm>

TECHNICAL SUPPORT

For technical support, please contact us at support@olink.com

For research use only. Not for use in diagnostic procedures.

This product includes a license for non-commercial use of Olink products. Commercial users may require additional licenses. Please contact Olink Proteomics AB for details.

There are no warranties, expressed or implied, which extend beyond this description. Olink

Proteomics AB is not liable for property damage, personal injury, or economic loss caused by this product.

Olink® is a registered trademark of Olink Proteomics AB.

Patents: See www.olink.com/patents

This product is sold under license from PHRI Properties, Inc. and may be used under PHRI Properties patent rights outside the field of human in vitro diagnostics.

Components in the Olink Probe Kit utilise Lightning-Link™ technology and are provided under license from Innova Biosciences.

© 2016-2019 Olink Proteomics AB. All third party trademarks are the property of their respective owners.