



Olink® Reveal Automation Guide



Day 1



Day 2



Reagent fill volumes		
Art. Nr.	Component	Volume (μL)
89001	Olink® Reveal PCR Solution	1520
89003	Olink® Reveal PCR Additive	1520
89006	Olink® Reveal Enzyme A	320
89007	Olink® Reveal Enzyme B	50
89008	Olink® Reveal Sample Buffer	1375
89009	Olink® Explore/Reveal Negative control	160
89010	Olink® Explore/Reveal Plate Control	160
89011	Olink® Explore/Reveal Sample Control	95

For full instructions on how to perform the Olink® Reveal workflow, see Olink® Reveal User Manual.