



Uppsala November 20th, 2025

Release Notes

NPX™ Map pre-processing software, ngs2counts, version 6.2.0

Release version v.6.2.0

Release Date: November 20th, 2025

Note: The ngs2counts v 6.2.0 software only supports pre-processing of data from Illumina sequencing instruments.

Description ngs2counts, version 6.2.0

This release introduces improved counts detection for Olink internal control assays, with enhancements relevant for **Olink® Explore HT** and **Olink® Reveal**.

There are no changes for Olink® Explore 384/3072.

Updating may affect NPX values, so changes should only be made between projects to maintain data consistency.

Note: For detailed information and guidance, please contact support@olink.com.

Content in ngs2counts version 6.2.0

Table 1 includes changes for ngs2counts version 6.2.0

Feature	Description and recommended action
Improved counts detection	Improved counts detection for Olink internal control assays, with enhancements relevant for Olink® Explore HT and Olink® Reveal . There are no changes for Olink® Explore 384/3072.

Software download

The new version of ngs2counts is available to download from the Olink® Software webpage:

<https://olink.com/software/download> Password: bg7Y0??yg89&!



Customer documentation

The Technical information NPX™ Map pre-processing, version 6.2.0 is available with corresponding changes from <https://olink.com/software/download> password: bg7Y0??yg89&! and on Document download center at Olink.com <https://www.olink.com/resources-support/document-download-center/>