



Olink[®] Flex — Flexibility without compromise

Olink Flex is a customizable made-to-order product that enables you to create tailored solutions to fit your unique inflammation research needs. Freely mix and match 5-30 biomarkers into a single panel with no compromise on data quality.

Mix and match

Design your own unique protein biomarker panel using Olink Insight

5-30

Combine 5-30 protein biomarkers and run 40 samples per kit

~200

Choose from a library of close to 200 pre-validated protein biomarkers

99%

Build your Flex panel with true 99% combinability of protein biomarkers

pg/mL

Readout in both absolute and relative quantification

Validation Data

Fully validated protein biomarker library

Build your Olink Flex panel from pre-validated protein biomarker assays using our online panel design tool in Olink[®] Insight. Each protein biomarker assay maintains high quality, performance, specificity and sensitivity for robust results.

To learn more visit us at olink.com or scan the QR code.

www.olink.com

© 2024 Olink Proteomics AB, part of Thermo Fisher Scientific.
Olink products and services are For Research Use Only. Not for use in diagnostic procedures.
This document is not intended to convey any warranties or representations of any kind.
1326, v1.2, 2024-09-30

